

Anti-Phospho-NF-kB p65 (Thr254) Rabbit antibody



Catalog Number: 166881

Product name

Anti-Phospho-NF-kB p65 (Thr254) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospho-NF-kB p65 (Thr254) Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around phosphorylation site of threonine 254 (F-R-T(p)-P-P) derived from Human NFκB-p65.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF


Dilutions

WB:1/500-1/1000 ICC:1/100-1/200 IHC:1/50-1/100


Validations



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using NFκB-p65 (Phospho-Thr254) Antibody (left) or the same antibody preincubated with blocking peptide (right).



Western blot analysis of extracts from 3T3 cells untreated or treated with TNF-α using NFκB-p65 (Phospho-Thr254) Antibody.




Western blot detection of NFκB-p65(Phospho-Ser254) in 3T3 cells untreated or treated with Alkaline phosphatase using NFκB-p65(Phospho-Ser254) Rabbit pAb (dilution 1:500, upper) or β-Actin Mouse mAb (200068-8F10, lower). Predicted band size:65kDa.Observed band size:65KDa.

Anti-Phospho-NF-kB p65 (Thr254) Rabbit antibody



Catalog Number: 166881



Immunofluorescence staining of methanol-fixed HeLa using NFκB-p65(Phospho-Thr254) Antibody .