

Human PPIL1/CYPL1/CGI-124/UNQ2425 (His tag) recombinant protein



Catalog Number: 518881

General Information

Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase-Like 1,
PPIase, Rotamase PPIL1, PPIL1, CYPL1

Protein Construction

A DNA sequence encoding the Human
PPIL1/CYPL1/CGI-124/UNQ2425 Met1-Gly166 is
expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
TrisHCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PPIL1/CYPL1/CGI-124/UNQ2425 (His tag)
recombinant protein