

# DBT antibody



Catalog Number: 109881

## Product name

DBT antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DBT Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A375 cells, K-562 cells. Observed molecular weight by Western-blot: 53kd

## Preparation

This antibody was obtained by immunization of DBT recombinant protein (Accession Number: NM\_001918). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109881(DBT antibody) at dilution of 1:300