

Nicastrin antibody



Catalog Number: 113181

Product name

Nicastrin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Nicastrin Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, C6 cells, human brain tissue. Positive IF detected in C6 cells. Observed molecular weight by Western-blot: 110kd

Preparation

This antibody was obtained by immunization of Nicastrin recombinant protein (Accession Number: NM_001290184). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113181(NCSTN antibody) at dilution of 1:800



Immunofluorescent analysis of C6 cells using Catalog No:113181(Nicastrin Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)