

Anti-FANCC antibody



Catalog Number: 105181

Product name

Anti-FANCC antibody

IHC-P: 0.1-2 µg/mL

ICC/IF: 0.5-1.5 µg/mL


Specificity

Human FANCC

Validations

Antibody description

Rabbit polyclonal to FANCC



Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human FANCC, and purified by antigen affinity chromatography.

FANCC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Formulation

0.2 µm filtered solution in PBS


Immunochemical staining of human FANCC in human brain with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.



sections).

FANCC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunochemical staining of human FANCC in human pancreas with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).

Applications


IHC-P, IF, ICC/IF

Dilutions

Anti-FANCC antibody



Catalog Number: 105181



FANCC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence

Immunofluorescence staining of FANCC in HeLa cells. Cells were fixed with 4% PFA, blocked with 10% serum. Then incubated with rabbit anti-human FANCC polyclonal antibody (1 $\mu\text{g/ml}$) (Figure A), incubated with rabbit anti-human FANCC polyclonal antibody and antigen (Figure B) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to cytoplasm and nucleus.