

# Anti-IGFBP6 Rabbit antibody



Catalog Number: 167181

## Product name

Anti-IGFBP6 Rabbit antibody

## Specificity

Human

## Antibody description

IGFBP6 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human IGFBP6

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/1000

ICC/IF: 1/20

FC: 1/100

## Validations



Western blot detection of IGFBP6 in K562, HeLa cell lysates using IGFBP6 Rabbit pAb(1:1000 diluted). Predicted band size: 25kDa. Observed band size: 35kDa.