

FAM186B antibody



Catalog Number: 110581

Product name

FAM186B antibody

Specificity

Human; other species not tested.

Antibody description

FAM186B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 120 kd

Preparation

This antibody was obtained by immunization of FAM186B recombinant protein (Accession Number: XM_006719625). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110581(FAM186B Antibody) at dilution of 1:600