

IAH1 antibody

Catalog Number: 111581

Product name

IAH1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

IAH1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue, mouse kidney tissue. Positive IP detected in mouse kidney tissue, mouse testis tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of IAH1 recombinant protein (Accession Number: NM_001039613). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111581(IAH1 antibody) at dilution of 1:500



IP Result of anti-IAH1 (IP:Catalog No:111581, 3ug; Detection:Catalog No:111581 1:500) with mouse kidney tissue lysate 5000ug.