

PAR2 antibody

Catalog Number: 113581

Product name

PAR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PAR2 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue. Observed molecular weight by Western-blot: 50-55kd

Preparation

This antibody was obtained by immunization of PAR2 recombinant protein (Accession Number: NM_005242). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations

mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:113581(F2RL1 antibody) at dilution of 1:300