

# CREBBP Polyclonal Antibody



Catalog Number: 161981

---

**Product name**

CREBBP Polyclonal Antibody

**Specificity**

Mouse, Rat

**Antibody description**

Polyclonal antibody to CREBBP

**Preparation**

Antigen: Recombinant protein of mouse CREBBP

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

WB, IHC

**Dilutions**

WB 1:500 - 1:2000

IHC 1:50 - 1:200

**Validations**

Immunohistochemistry - CREBBP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using CREBBP antibody at dilution of 1:100 (40x lens).