

PNPLA3 antibody

Catalog Number: 113981

Product name

PNPLA3 antibody

Specificity

Human; other species not tested.

Antibody description

PNPLA3 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, human adipose tissue, Jurkat cells. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of PNPLA3 recombinant protein (Accession Number: BC014449). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113981(PNPLA3 antibody) at dilution of 1:100



IP Result of anti-PNPLA3 (IP:Catalog No:113981, 4ug; Detection:Catalog No:113981 1:100) with HepG2 cells lysate 1000ug.