

# TMPRSS11A Polyclonal Antibody



Catalog Number: 165981

## Product name

Rabbit IgG

TMPRSS11A Polyclonal Antibody

## Applications

## Specificity

WB

Human, Mouse, Rat

## Dilutions

## Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to TMPRSS11A

## Validations

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 300-421 of human TMPRSS11A (NP\_872412.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type



## Western blot - TMPRSS11A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TMPRSS11A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.