

Bersanlimab reference antibody



Catalog Number: 176981

Product name

Bersanlimab reference antibody

Specificity

Human ICAM1 / CD54

Antibody description

Anti-ICAM1 / CD54 Antibody (Bersanlimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SEC-HPLC

The purity of Anti-ICAM1 / CD54 Antibody (bersanlimab) is 100%, determined by SEC-HPLC.



ELISA

Immobilized human ICAM-1 / CD54 Protein, His Tag Protein at 2 ug/mL can bind anti-ICAM1 / CD54 Antibody (bersanlimab) $EC_{50}=0.002577$



ug/mL.

SDS-PAGE

Anti-ICAM1 / CD54 Antibody (bersanlimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95.0%.