

hCG alpha antibody

Catalog Number: 111281

Product name

hCG alpha antibody

Specificity

Human; other species not tested.

Antibody description

hCG alpha Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human placenta tissue, human lung cancer tissue.

Preparation

This antibody was obtained by immunization of hCG alpha recombinant protein (Accession Number: NM_000735). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111281(CGA Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:111281(CGA Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:111281(CGA Antibody) at dilution of 1:200 (under 40x lens).