

Human TMEFF1/Tomoregulin-1 (Fc Tag) recombinant protein



Catalog Number: 502281

General Information

Gene Name Synonym

H7365; Transmembrane protein with EGF-like and one follistatin-like domain

Protein Construction

A DNA sequence encoding the human TMEFF1 (NP_003683.2) (Met1-Val330) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 40

Molecular Mass

The recombinant human TMEFF1 consists 532 amino acids and predicts a molecular mass of 58.5 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human TMEFF1 / Tomoregulin-1 Protein (Fc Tag)