

# TMTC1 antibody



Catalog Number: 116281

## Product name

TMTC1 antibody

## Immunogen

[Human TMTC1 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse; other species not tested.

## Antibody description

TMTC1 Rabbit Polyclonal antibody. Positive WB detected in mouse bladder tissue, mouse spleen tissue, mouse colon tissue. Observed molecular weight by Western-blot: 99 kDa, 65 kDa

## Preparation

This antibody was obtained by immunization of TMTC1 recombinant protein (Accession Number: NM\_175861). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse bladder tissue were subjected to SDS PAGE followed by western blot with Catalog No:116281(TMTC1 antibody) at dilution of 1:1000