

Anti-NOX4 antibody



Catalog Number: 142681

Product name

Anti-NOX4 antibody

Specificity

Human NOX4

Antibody description

Rabbit Polyclonal to Human NOX4

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human NOX4, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 1-4 µl/mg of lysate

Validations



Human NOX4 Immunoprecipitation(IP) 14669