

EIF5B Polyclonal Antibody



Catalog Number: 164681

Product name

EIF5B Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to EIF5B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 901-1220 of human EIF5B (NP_056988.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:10 - 1:100

Validations



Immunofluorescence - EIF5B Polyclonal Antibody

Immunofluorescence analysis of A549 cells using EIF5B antibody . Blue: DAPI for nuclear staining.