

# Anti-Peroxiredoxin 1/PAG Rabbit antibody



Catalog Number: 167681

## Product name

Rabbit IgG

Anti-Peroxiredoxin 1/PAG Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/2000-1/10000

Peroxiredoxin 1/PAG Rabbit monoclonal antibody

ICC/IF: 1/20

## Preparation

FC: 1/20

Antigen: A synthetic peptide of human  
Peroxiredoxin 1/PAG

IP: 1/20

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage


Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of Peroxiredoxin 1/PAG in  
A549,HL-60 using Peroxiredoxin 1/PAG Rabbit  
mAb(1:1000 diluted)