

ABCA7 antibody

Catalog Number: 107681

Product name

ABCA7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ABCA7 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 200-240 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_019112). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107681(ABCA7 Antibody) at dilution of 1:1000