

# NUP62 antibody

Catalog Number: 113381

## Product name

NUP62 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NUP62 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, NIH/3T3 cells. Positive IP detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 62kd

## Preparation

This antibody was obtained by immunization of NUP62 recombinant protein (Accession Number: NM\_153719). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113381(NUP62 antibody) at dilution of 1:1500



IP Result of anti-NUP62 (IP:Catalog No:113381, 4ug; Detection:Catalog No:113381 1:1000) with NIH/3T3 cells lysate 1200ug.