

PSMB11 antibody



Catalog Number: 114381

Product name

PSMB11 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PSMB11 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 32-35kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001099780). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114381(PSMB11 antibody) at dilution of 1:300