

# Anti-SLC39A9 antibody



Catalog Number: 176381

## Product name

Anti-SLC39A9 antibody

## Specificity

Human, Rat

## Antibody description

Rabbit polyclonal antibody to SLC39A9

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human slc39a9. aa range:101-150

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations



Fig1.; Western Blot analysis of various cells using ZIP9 Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000