

IL-16 antibody



Catalog Number: 107381

Product name

IL-16 antibody

Specificity

Human; other species not tested.

Antibody description

IL-16 Mouse Monoclonal antibody. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of IL-16 recombinant protein (Accession Number: BC040272). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107381(IL16 antibody) at dilution of 1:20000