

# Collagen Type IV antibody



Catalog Number: 109381

## Product name

Collagen Type IV antibody

IP: 1:200-1:1000

IF: 1:10-1:100


## Specificity

Human; other species not tested.

## Validations

## Antibody description

Collagen Type IV Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in HeLa cells, A431 cells, SGC-7901 cells. Observed molecular weight by Western-blot: 50~55 kDa




## Preparation

This antibody was obtained by immunization of Collagen Type IV recombinant protein (Accession Number: BC005305). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109381(COL4A6 antibody) at dilution of 1:1000

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

IP Result of anti-COL4A6 (IP:Catalog No:109381, 4ug; Detection:Catalog No:109381 1:300) with mouse skeletal muscle tissue lysate 3600ug.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IF, IP

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of HeLa cells using Catalog No:109381(COL4A6 Antibody) at dilution

# Collagen Type IV antibody



Catalog Number: 109381

---

of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)