

## General Information

### Gene Name Synonym

Thy-1 membrane glycoprotein variant 2; cDNA FLJ33325 fis, clone BNGH42008832, highly similar to THY-1 MEMBRANE GLYCOPROTEIN; cDNA FLJ33677 fis, clone BRAWH2002125, highly similar to THY-1 MEMBRANE GLYCOPROTEIN

### Protein Construction

A DNA sequence encoding the human THY1 (NP\_006279.2) (Met1-Cys130) was expressed with five amino acids (DDDDK) at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per µg protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Gln 20

## Molecular Mass

The recombinant human THY1 consists of 117 amino acids and predicts a molecular mass of 13.2 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

### SDS-PAGE



Human CD90/THY-1 Protein 15092