

# Anti-Adipose Triglyceride Lipase Rabbit antibody



Catalog Number: 166781

## Product name

WB: 1/1000

Anti-Adipose Triglyceride Lipase Rabbit antibody

IHC: 1/500

## Specificity

ICC/IF: 1/100

Human, Mouse, Rat

IP: 1/20

## Validations

## Antibody description

Adipose Triglyceride Lipase Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Adipose Triglyceride Lipase

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions



Immunohistochemistry of Adipose Triglyceride Lipase in paraffin-embedded Human breast cancer tissue using Adipose Triglyceride Lipase Rabbit pAb at dilution 1/100



Western blot detection of Adipose Triglyceride Lipase in K562, C6, 3T3, HeLa cell lysates using Adipose Triglyceride Lipase Rabbit pAb(1:1000 diluted). Predicted band size: 55kDa. Observed band size: 55kDa.