

# Anti-BAG2 Rabbit antibody



Catalog Number: 168781

## Product name

Anti-BAG2 Rabbit antibody

## Specificity

Human

## Antibody description

BAG2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human BAG2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB, FC, IP

## Dilutions

WB: 1/1000

## Validations



Western blot detection of BAG2 in Jurkat, HeLa cell lysates using BAG2 Rabbit mAb(1:1000 diluted). Predicted band size: 24kDa. Observed band size: 24kDa.