

C6orf173 antibody



Catalog Number: 108781

Product name

C6orf173 antibody

Immunogen

[Human C6orf173 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C6orf173 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, A549 cells. Observed molecular weight by Western-blot: 7-10 kDa

Preparation

This antibody was obtained by immunization of C6orf173 recombinant protein (Accession Number: NM_001012507). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:108781(C6orf173 Antibody) at dilution of 1:600