

# H3N2 HA (His Tag) recombinant protein



Catalog Number: 500091

## General Information

### Gene Name Synonym

Hemagglutinin HA1 chain; Hemagglutinin HA2 chain

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/Hong Kong/1/1968(H3N2)) hemagglutinin (Q91MA7) (Met1-Trp530) was expressed with a C-terminal polyhistidine tag.

### Organism

H3N2

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

1. Measured by its ability to bind with Neu5Aca2-3Galb1-4GlcNAcb-PAA-biotin (01-077) using the Octet RED System.
2. Measured by its ability to bind with Neu5Aca2-6GalNAca-PAA-biotin (01-059) using the Octet RED System.
3. Measured by its ability to agglutinate guinea pig red blood cells. HA titer is 0.1-0.6 µg/mL for 1% GRBC.

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

## Predicted N terminal

Gln 17

## Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Hong Kong/1/1968(H3N2)) comprises 525 amino acids and has a predicted molecular mass of 59 kDa. The apparent molecular mass of the protein is approximately 65 kDa in SDS-PAGE under reducing conditions.

## Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

## SDS-PAGE



# H3N2 HA (His Tag) recombinant protein



Catalog Number: 500091

---

Influenza A H3N2 (A/Hong Kong/1/1968)

Hemagglutinin / HA Protein (His Tag) SDS-PAGE