

PGRMC2 Polyclonal Antibody



Catalog Number: 163091

Product name

PGRMC2 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PGRMC2

Preparation

Antigen: A synthetic peptide of human PGRMC2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IF

Dilutions


WB 1:500 - 1:2000
IF 1:50 - 1:200

Validations



Western blot - PGRMC2 Polyclonal Antibody

Western blot analysis of extracts of MCF-7 cells, using PGRMC2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5min.



Immunofluorescence - PGRMC2 Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using PGRMC2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.