

PPP1R7 antibody



Catalog Number: 114091

Product name

PPP1R7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PPP1R7 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HEK-293 cells, PC-3 cells. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of PPP1R7 recombinant protein (Accession Number: NM_002712). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114091(PPP1R7 Antibody) at dilution of 1:600