

Anti-Phospho-ERK1/2(Tyr222/205) Rabbit antibody



Catalog Number: 169091

Product name

Polyclonal

Anti-Phospho-ERK1/2(Tyr222/205) Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

ICC/IF, WB, IHC-P, ELISA

Antibody description

Phospho-ERK1/2(Tyr222/205) Rabbit polyclonal antibody

Dilutions

IF: 1/50-200

Preparation

Antigen: Synthesized phospho-peptide around the phosphorylation site of human ERK 1/2 (phospho Tyr222/205)

WB: 1/500-2000

IHC: 1/50-300

IHC: 1/50-300

Formulation


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Validations



Western Blot analysis of various cells using Phospho-ERK 1/2 (Y222/205) Polyclonal Antibody