

EVI1 antibody

Catalog Number: 110491

Product name

EVI1 antibody

Specificity

Human; other species not tested.

Antibody description

EVI1 Rabbit Polyclonal antibody. Positive IF detected in SKOV-3 cells. Positive WB detected in SKOV-3 cells. Observed molecular weight by Western-blot: 118 kDa

Preparation

This antibody was obtained by immunization of EVI1 recombinant protein (Accession Number: XM_005247220). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:50-1:500

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110491(EVI1 Antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed SKOV-3 cells using Catalog No:110491(EVI1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)