

Anti-DLEU7 antibody



Catalog Number: 176491

Product name

Polyclonal

Anti-DLEU7 antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

IHC-P

Rabbit polyclonal antibody to DLEU7

Dilutions

Preparation

IHC-P: 1:50-200

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human dleu7. aa range:51-100

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1:: Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200