

TULP1 antibody

Catalog Number: 116491

Product name

TULP1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

human, mouse, rat; other species not tested.

Validations

Antibody description

TULP1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IF detected in Hela. Positive IHC detected in human eye. Observed molecular weight by Western-blot: 70kd




Preparation

This antibody was obtained by immunization of TULP1 recombinant protein (Accession Number: NM_001289395). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116491(TULP1 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human eye using Catalog No:116491(TULP1 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human eye using Catalog No:116491(TULP1

TULP1 antibody

Catalog Number: 116491

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using TULP1 antibody Catalog No:116491 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).