

Rab18 antibody



Catalog Number: 107491

Product name

Rab18 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Rab18 Mouse Monoclonal antibody. Positive WB detected in Human adrenal gland, mouse testis tissue. Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of Rab18 recombinant protein (Accession Number: NM_021252). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:107491(RAB18 antibody) at dilution of 1:1000