

Mouse Cndp2/Cn2 (His tag) recombinant protein



Catalog Number: 519491

General Information

Gene Name Synonym

CNDP2/CNDP dipeptidase 2/Cytosolic non-specific dipeptidase/Glutamate carboxypeptidase-like protein 1/Cn2

Protein Construction

A DNA sequence encoding the Mouse Cndp2/Cn2 Met1-Asn475 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Cndp2/Cn2 (His tag) recombinant protein