

JNK antibody

Catalog Number: 111891

Product name

JNK antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:500-1:5000


IHC: 1:20-1:200

Antibody description

JNK Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human lung cancer tissue, mouse heart tissue. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, A431 cells, HEK-293 cells, Jurkat cells, PC-3 cells, rat spleen tissue, RAW264.7 cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 40-45 kDa, 50-55 kDa

IF: 1:20-1:200

Validations



Preparation

This antibody was obtained by immunization of JNK recombinant protein (Accession Number: NM_001323323). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111891(JNK antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF



Dilutions



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:111891(JNK antibody) at dilution of 1:50 (under 10x lens)

JNK antibody

Catalog Number: 111891



IP Result of anti-JNK (IP:Catalog No:111891, 3ug; Detection:Catalog No:111891 1:1000) with HeLa cells lysate 1000ug.

Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:111891(JNK Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)