

AHSP Polyclonal Antibody



Catalog Number: 164891

Product name

AHSP Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to AHSP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human AHSP (NP_057717.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Western blot - AHSP Polyclonal Antibody