

Anti-LAMP2 Rabbit antibody



Catalog Number: 167891

Product name

Anti-LAMP2 Rabbit antibody

Specificity

Human

Antibody description

LAMP2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human LAMP2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/200

Validations



Western blot detection of LAMP2 in HeLa,A549,U2OS,U87-MG cell lysates using LAMP2 Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:110kDa.