

Human ACADM (His tag) recombinant protein



Catalog Number: 519891

General Information

Gene Name Synonym

Medium-Chain Specific Acyl-CoA Dehydrogenase
Mitochondrial, MCAD, ACADM

Protein Construction

A DNA sequence encoding the Human ACADM
Lys26-Asn421 is expressed with a 6His tag at the
N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM
Tris, 0.1M NaCl, 20% Glycerol, pH 8.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ACADM (His tag) recombinant protein