

EIF3C antibody

Catalog Number: 110191

Product name

EIF3C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EIF3C Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, MCF7 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of EIF3C recombinant protein (Accession Number: NM_003752). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110191(EIF3C antibody) at dilution of 1:500



IP Result of anti-EIF3C (IP:Catalog No:110191, 4ug; Detection:Catalog No:110191 1:500) with K-562 cells lysate 3200ug.