

NISCH antibody

Catalog Number: 113191

Product name

NISCH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NISCH Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 190-200 kDa

Preparation

This antibody was obtained by immunization of NISCH recombinant protein (Accession Number: NM_007184). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113191(NISCH antibody) at dilution of 1:300



IP Result of anti-NISCH (IP:Catalog No:113191, 4ug; Detection:Catalog No:113191 1:500) with mouse brain tissue lysate 4000ug.