

# ACTL6B Polyclonal Antibody



Catalog Number: 165191

## Product name

ACTL6B Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ACTL6B

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ACTL6B (NP\_057272.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

## Validations




Western blot - ACTL6B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ACTL6B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 40s.



Immunohistochemistry - ACTL6B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using ACTL6B antibody at dilution



of 1:100 (40x lens).

# ACTL6B Polyclonal Antibody




Catalog Number: 165191

---


Immunohistochemistry - ACTL6B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney cancer using ACTL6B antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ACTL6B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse jejunum using ACTL6B antibody at dilution of 1:100 (40x lens).



Immunofluorescence - ACTL6B Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using ACTL6B antibody .