

# TMEM54 antibody



Catalog Number: 116191

## Product name

TMEM54 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TMEM54 Rabbit Polyclonal antibody. Positive WB detected in SKOV-3 cells. Observed molecular weight by Western-blot: 24 kDa; 43 kDa

## Preparation

This antibody was obtained by immunization of TMEM54 recombinant protein (Accession Number: NM\_033504). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116191(TMEM54 Antibody) at dilution of 1:300