

Human PAFAH1B2/PAFAHB (His tag) recombinant protein



Catalog Number: 519191

General Information

Gene Name Synonym

Platelet-Activating Factor Acetylhydrolase IB
Subunit Beta, PAF Acetylhydrolase 30 kDa
Subunit, PAF-AH 30 kDa Subunit, PAF-AH Subunit
Beta, PAFAH Subunit Beta, PAFAH1B2, PAFAHB

Protein Construction

A DNA sequence encoding the Human
PAFAH1B2/PAFAHB Ser2-Ala229 is expressed
with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM
Tris, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PAFAH1B2/PAFAHB (His tag)
recombinant protein