

# FAM48A antibody



Catalog Number: 110591

## Product name

FAM48A antibody

## Specificity

Human; other species not tested.

## Antibody description

FAM48A Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, Jurkat cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 86 kDa

## Preparation

This antibody was obtained by immunization of FAM48A recombinant protein (Accession Number: NM\_001014286). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110591(FAM48A Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110591(FAM48A Antibody) at dilution of 1:50 (under 40x lens)