

# RBBP5 antibody



Catalog Number: 114591

---

## Product name

RBBP5 antibody

## Immunogen

[Human RBBP5 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RBBP5 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 60-65 kDa

## Preparation

This antibody was obtained by immunization of RBBP5 recombinant protein (Accession Number: NM\_001193272). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114591(RBBP5 antibody) at dilution of 1:400