

Anti-Ly-6A/E Rabbit antibody



Catalog Number: 167591

Product name

Anti-Ly-6A/E Rabbit antibody

Specificity

Mouse

Antibody description

Ly-6A/E Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of mouse Ly-6A/E

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/50

Validations



Western blot detection of Ly-6A/E in HeLa,A549,MCF-7,HL-60 cell lysates using Ly-6A/E Rabbit mAb(1:1000 diluted).Predicted band size:14kDa.Observed band size:18kDa.